宁波新跃医疗科技股份有限公司
NINGBO XINWELL MEDICAL TECHNOLOGY CO., LTD.
 使用说明书
 Instruction for Use
 Gebruiksaanwijzing

产品名称: 一次性使用结扎夹
 Product name Disposable ligation clips
 Productnaam Wegwerpliggingsklem


文件编号: RD/S-105-0002
 Doc. No.
 Doc. Nee

文件版本: A/0
 Doc. Version
 Doc. Version

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3.0	No.188 Binjiang Road,Cixi High-tech Industrial Development Zone, Cixi City, Zhejiang, 315301, China	任振俦 Zhenchu Ren	2020-12-11	DHF-105-40-07
4.0	The clip body should have a complete shape and a smooth feel, except for the part where the inner side of the product's arc nail leg is in contact with human tissue. The product should be free of defects such as burrs, flashes, cracks, and scars. Het omhulsel van deze clip moet volledig glad aanvoelen, met uitzondering van het deel aan de binnenkant van de boog nagelpoot dat in contact komt met menselijk weefsel. Het product mag geen bramen, scheuren en andere schade bevatten en moet ook vrij zijn van defecten.	任振俦 Zhenchu Ren	2021-01-20	ECN2021012001
A/0	Change version 4.0 to A/0	任振俦 Zhenchu Ren	2021-04-02	ECN2021040202

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[Parts Description]

The disposable ligation clip has a two-arm arc structure having a "V" shape, which is composed of a locking mechanism (left and right), hinge, bow-shape legs, integrated ridges and bosses (left and right). The structure diagram is shown in Fig. 1.


Fig. 1

[Product Description]

The disposable ligation clip is a non-absorbable, non-active implantable and single use medical device designed for use in general surgical procedures that require vessel or bile duct ligation. The clip is made of medical grade Homopolyformaldehyde (POM).

[Performance Index]

1. Appearance

1) The clip body should have a complete shape and a smooth feel, except for the part where the inner side of the product's arc nail leg is in contact with human tissue. The product should be free of defects such as burrs, flashes, cracks, and scars.

2) The ridges shall be on the inner side of the bow shape legs. Ridges shall be clear and complete, free of defects such as overflow angle and sharp edges.

2. Physical performance

1) Reliability of the lock: The lock shall have good reliability. There shall be no popping or slipping after the lock.

2) Toughness: The ligation clip should have good toughness. After clamping, no cracks or fractures should occur.

3) Clamping performance: The device shall fit well with the clip applier and have good clamping performance.

3. Chemical performance

Ethylene oxide residue: The ethylene oxide residue of the device shall be not more than 10ug/cm².

4. Biological performance

Sterility: The device shall be sterile after sterilization with ethylene oxide.

[Intended Use]

Disposable ligation clip is intended for use in general surgical procedures involving permanent ligation of vessels or bile duct.

[Contraindication]

- Disposable ligation clip is not intended for use as a fallopian contraceptive tubal occlusion device.
- Disposable ligation clip is contraindicated for use in ligating the renal artery during laparoscopic donor nephrectomies.

Attention

1. Devices used in combination with disposable ligation clip

1) Clip Applier

- The clip applier is reusable which manually load and deliver one clip at a time.
- Alignment of the jaw is critical for safe application of the clip. Always check the alignment of the applier jaws before use.
- It is recommended to use "Xinwell" brand clip appliers. Be sure to verify that the specifications of the clip appliers used to match with the specified ligation clip!

2) Clip Remover (if applicable)

- The clip remover is reusable which manually remove the clip at one time.
- Alignment of the jaw is critical for safe application of the clip. Always check the alignment of the clip remover jaws before use.
- It is recommended to use "Hem-o-lok" brand clip removers. Be sure to verify that the specifications of the clip removers used to match with the specified ligation clip!

2. Before applying a clip, verify the structural size and condition of the vessel or structure and use the proper size clip.

3. Before using the ligation clip, make sure that the ligation clip is in place.

4. Security of the closure should be confirmed after ligation.
5. Application of a second clip on all other vessels should be dictated by surgeons.

 Warning

1. Do not use if package is damaged!
2. For single use only. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure.
3. Do not use if it exceeds the expiry date.
4. It shall be destroyed and disposed in accordance with the relevant laws and regulations of relevant countries, and it is strictly prohibited to discard it at will.
5. Sterile pack should be used up at one time. Please destroy the unused ones!
6. This product must be operated by a professional surgeon.

 MR Safe

Disposable ligation clips are constructed of a non-absorbable polymer material and are considered MR Safe according to the current criteria and labeling terminology, (MR safe is an item that poses no known hazards in all MRI environments. Using the terminology, “MR Safe” items are non-conducting, nonmetallic and non-magnetic items).

[Usage]

● Loading clips

1. To load the applier, grasp the applier and carefully insert the jaws of the applier into the cartridge slot, making sure the jaws are perpendicular to the base of the cartridge. Gently press the applier over the clip until there is an audible double click (Fig. 2). Do not force the applier into the cartridge or onto the clip. The applier should enter and withdraw from the cartridge easily.
2. Remove the applier from the cartridge ensuring the clip is held securely in the applier jaws. The clip bosses should seat in the notches of the applier jaws (Fig. 3). It may be necessary to hold the cartridge to allow the clip to be removed.
3. Compress the applier handles and choose an appropriate way (open surgery or through Trocar) to enter the body.

● Clip positioning and closure

4. Orient the single tooth of the clip. This allows the user to visually confirm encapsulation of the structure being ligated. Note: ligate at a distance of 2-3mm from the distal of the tissue (or at least 4-6mm between the two ligature clips) in order to leave a margin. (Fig 4).
5. Apply sufficient force to the applier handles so the jaws close and the clip locks shut (Fig 5). Releasing pressure on the applier handles allows the applier to return to a fully open position (Fig 6).


Fig. 2


Fig. 3

Fig. 4


Fig. 5


Fig. 6

● Removing clips (if applicable)

6. Apply sufficient force to the applier handles so the jaws close and the clip locks shut (Fig 5). Releasing pressure on the applier handles allows the applier to return to a fully open position (Fig 7).

7. Insert the clip remover to enter the body . The clip shall be approached from the hinge side, not the locking mechanism side. Make sure that the jaws are directly over and in line with the legs of the clip (Fig. 7). Advance the remover until the user visualizes the hinge of the clip resting at the back of the remover jaws. NOTE: The hinge of the clip shall rest in this channel to ensure disengagement of the clip legs.

8. Once the hinge is located in this channel, apply steady pressure to the remover handles until a slight click is felt (Fig. 8). This indicates the clip legs have disengaged (Fig. 9). Visually confirm the clip has opened sufficiently to allow removal.


Fig 7


Fig 8


Fig 9


[Transportation and storage conditions]

The product shall be stored in a clean room that is light proof, rat proof, fireproof, insect proof, well ventilated and free of corrosive gas.

Ambient Temperature: - 5 °C ~ 40 °C

Relative Humidity: 10% ~ 80%

[Packaging marks]

No.	Symbol	Meaning	No.	Symbol	Meaning
1		DATE OF MANUFACTURE	8		STERILIZED USING ETHYLENE OXIDE
2		DO NOT REUSE	9		CONSULT INSTRUCTIONS FOR USE
3		CAUTION	10		USE-BY DATE
4		MANUFACTURER	11		BATCH CODE
5		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY	12		DO NOT USE IF PACKAGE IS DAMAGED
6		THE PRODUCT HAS PASSED CE CERTIFICATION	13		CATALOGUE NUMBER
7		DO NOT RESTERILIZE	14		KEEP AWAY FROM SUNLIGHT
			15		KEEP DRY

The device is sterilized by ethylene oxide, and the shelf life is 3 years after sterilization.


Ningbo Xinwell Medical Technology Co., Ltd

Add: No.188 Binjiang Road,Cixi High-tech Industrial Development Zone, Cixi City, Zhejiang, 315301, China

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