

Disposable ligation clips



New Technology For Well Being

Products Features

Safety Locking Machanism

360° all-round ligation, and the lock is safe and reliable; the hand feels clear when closed, and the safety information is transmitted.

Bow-shape legs


The arc design has a larger ligation range; The force is even and there is no cutting effect.

2 Carina

Accurate clamping and not easy to fall off; accurate guidance helps efficient ligation.

4 Hinge

Adjust the hinge according to the force; if the ligation fails ,it will bounce off clearly.



5 Integrated ridges

The curved non-slip surface is more consistent with the blood vessel interface; The integrated non-slip tooth effectively prevents sliding.

	Model	The size range of closed blood vessels and cystic duct	Size
Manual arx	JZJ-ML-1	3~10mm	Medium Large (ML)
	JZJ-ML-2	3~10mm	Medium Large (ML)
	JZJ-ML-3	3~10mm	Medium Large (ML)
	JZJ-ML-4	3~10mm	Medium Large (ML)
	JZJ-ML-5	3~10mm	Medium Large (ML)
	JZJ-ML-6	3~10mm	Medium Large (ML)
Y In wall accor	JZJ-L-1	5~13mm	Large (L)
	JZJ-L-2	5~13mm	Large (L)
	JZJ-L-3	5~13mm	Large (L)
	JZJ-L-4	5~13mm	Large (L)
	JZJ-L-5	5~13mm	Large (L)
	JZJ-L-6	5~13mm	Large (L)
Xinuel inix	JZJ-XL-1	7~16mm	Extra Large (XL)
	JZJ-XL-2	7~16mm	Extra Large (XL)
	JZJ-XL-3	7~16mm	Extra Large (XL)
	JZJ-XL-4	7~16mm	Extra Large (XL)
	JZJ-XL-5	7~16mm	Extra Large (XL)
	JZJ-XL-6	7~16mm	Extra Large (XL)

Ningbo Xinwell Medical Technology Co., Ltd.

No.188 Binjiang Road, Cixi High-tech Industrial Development Zone, Cixi City, Zhejiang, P.R. China

Web: www.xinwell.com

P.C.: 315301

